

Rabbit Anti-PGAM1 antibody

SLM-60524R

Product Name	PGAM1
Chinese Name	
Immunogen Species	Rabbit
Clonality	Monoclonal
Clone NO.	D3C4
React Species	(predicted: Human, Rat,) WB=1:500-2000,IHC-P=1:50-200,IHC-F=1:50-200,IF=1:50-200,Flow-Cyt=1:50-100
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Cellular localization	The cell membrane Extracellular matrix Secretory protein
Form	Liquid
Concentration	1mg/ml
Lsotype	IgG/Kappa
Purification	Affinity purified by Protein A
Buffer Solution	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	以 2, 3-二磷酸甘油酯为反应起始剂的 3-和 2-磷酸甘油酯的相互转化。也能催化 EC5.4.2.4 (合酶) 的反应, 但活性降低。