

Mouse Anti-CEP55 antibody

SLM-60527M

Product Name	CEP55
Chinese Name	
Immunogen Species	Mouse
Clonality	Monoclonal
Clone NO.	B7F5
React Species	(predicted: Human,) WB=1:500-1000,ICC/IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Cellular localization	The cell membrane
Form	Liquid
Concentration	1mg/ml
Lsotype	IgG1/Kappa
Purification	Affinity purified by Protein G
Buffer Solution	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	在有丝分裂退出和胞质分裂中起作用。在胞质分裂期间，将 PDCD6IP 和 TSG101 招募到中体。成功完成胞质分裂所必需的。对大脑和肾脏的发育起作用。