



Mouse Anti-SNRPN antibody

SLM-60533M

Product Name	SNRPN
Chinese Name	
Immunogen Species	Mouse
Clonality	Monoclonal
Clone NO.	A3B1
React Species	(predicted: Human, Mouse, Rat,) WB=1:500-1000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Cellular localization	The nucleus
Form	Liquid
Concentration	1mg/ml
Lsotype	IgG1/Kappa
Purification	Affinity purified by Protein G
Buffer Solution	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	可能参与组织特异性替代 RNA 加工事件。