

Rabbit Anti-SOX11 antibody

SLM-60624R

| | |
|------------------------------|--|
| Product Name | SOX11 |
| Chinese Name | |
| Immunogen Species | Rabbit |
| Clonality | Monoclonal |
| Clone NO. | B10G9 |
| React Species | Human,Mouse,Rat |
| Applications | IHC-P=1:50-200,IHC-F=1:50-200,IF=1:50-200 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user. |
| Cellular localization | The nucleus |
| Form | Liquid |
| Concentration | 1mg/ml |
| Lsotype | IgG |
| Purification | Affinity purified by Protein A |
| Buffer Solution | 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Attention | This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. |
| PubMed | PubMed |
| Product Detail | 作为转录激活因子的转录因子。 |