

Rabbit Anti-FUBP1 antibody

SLM-60625R

Product Name FUBP1

Chinese Name

Immunogen Species Rabbit

Clonality Monoclonal

Clone NO. G5E4

React Species Human,Mouse,Rat

Applications WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:50-200,IF=1:100-500,Flow-Cy
(Paraffin sections need antigen repair)
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Cellular localization The nucleus

Form Liquid

Concentration 1mg/ml

Lsotype IgG

Purification Affinity purified by Protein A

Buffer Solution 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product Detail

远端上游原件 Binding protein1 (FUBP1)是一种单链 DNABinding protein,有时也能与单链 DNA 结合。FUBP1 是原癌基因 c-myc 转录过程的反式激活因子,即 FUBP1 结合到管家基因 c-myc 的启动子元件 FUSE 上,增强 c-myc 的转录,从而促进细胞的增殖,在 Tumour 形成中具有重要作用。