


Mouse Anti-Beta Actin (CAEEL ,Loading Control)antibody

SLM-8770M

Product Name	Beta Actin (CAEEL,Loading Control)
Chinese Name	线虫肌动蛋白（内参）单克隆抗体
Alias	act-1; Actin-1; ACT1_CAEEEL; C-actin; β -actin; β actin;
Product Type	Internal reference anti
Immunogen Species	Mouse
Clonality	Monoclonal
Clone NO.	T1F5
React Species	Caenorhabditis Elegans WB=1:500-2000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Theoretical molecular weight	42kDa
Cellular localization	cytoplasmic
Form	Liquid
Concentration	1mg/ml
immunogen	Recombinant Caenorhabditis elegans Actin.
Lsotype	IgG1
Purification	affinity purified by Protein G
Buffer Solution	Caenorhabditis Elegans1M TBS(pH7.4) with 1% BSA, Caenorhabditis Elegans3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	

Product Picture



Sample:

Lane 1: Caenorhabditis Elegans tissue Lysate

Primary: Anti-Beta Actin (CAEEL) (SLM-8770M) at 1/1000 dilution

Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution

Predicted band size: 42 kD

Observed band size: 42 kD