

Mouse Anti-Tubulin-alpha (CAEEL , Loading Control)antibody

SLM-8772M

Product Name	Tubulin-alpha (CAEEL, Loading Control)
Chinese Name	线虫 Tubulin-alpha (内参) 单克隆抗体
Alias	Tubulin alpha chain; tba-1; C-tubulin-alpha; α -tubulin;
Product Type	Internal reference anti
Immunogen Species	Mouse
Clonality	Monoclonal
Clone NO.	T1F7M
React Species	Caenorhabditis Elegans WB=1:500-2000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Theoretical molecular weight	51kDa
Cellular localization	cytoplasmic
Form	Liquid
Concentration	1mg/ml
immunogen	Recombinant Caenorhabditis elegans Tubulin-alpha.
Lsotype	IgG1
Purification	affinity purified by Protein G
Buffer Solution	Caenorhabditis Elegans1M TBS(pH7.4) with 1% BSA, Caenorhabditis Elegans3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	

Product Picture

