



Rabbit Anti-FITC / FITC antibody

bsF-0366R

Product Name	Rabbit Anti-FITC / FITC
Chinese Name	标记兔抗 FITC
Alias	fluorescein isothiocyanate; Fluorescein 5-isothiocyanate. 荧光素标记兔抗 FITC
Immunogen Species	Rabbit
Clonality	Polyclonal
Applications	IHC-P=1:20-200 IHC-F=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
Concentration immunogen	1.0 mg/ml FITC
Lsotype	
Buffer Solution	10 mM TBS (pH=7.4) with 1% BSA, 3% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
Product Detail	