

Eukaryotic Interleukin 1 Beta (IL1b)

Instruction Manual

SFPA147Hu61

Homo sapiens (Human)

Source	Eukaryotic expression
Host	293F cell
Endotoxin Level	<1.0EU per 1ug (determined by the LAL method)
Subcellular Location	Secreted, Lysosome, Cytoplasm
Predicted Molecular Mass	19.0kDa
Accurate Molecular Mass	19kDa(Analysis of differences refer to the manual)
Residues & Tags	Ala117~Ser269 with N-terminal His Tag
Buffer Formulation	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	5.9
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
APVR SLNCTLRDSQ QKSLVMSGPY ELKALHLQGQ  
DMEQQVVFMS SFVQGEESND KIPVALGLKE KNLYLSCVLK DDKPTLQLES  
VDPKNYPKKK MEKRFVFNKI EINNKLEFES AQFPNWIYST SQAENMPVFL  
GGTKGGQDIT DFTMQFVSS
```

USAGE

Reconstitute in PBS or others.

STORAGE


Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image

TCGACCTGTCGGATCCAGTCGACCTGGGGACTCACGCAAAAGCTTGGTGAATGCTGGTCCATGAGATCGAAGCTCTGGACCTGGGGACGGGATCGACGAGCTGGTCTCCGATGTCCTTGGACAGGAGAGCAAGTATGCAAAATACCTGGTGGCTTGGGCTCAGGCAAGAT
LFFFRSLRSLRSLRSLRSD50QKSLVMSGPFYELKALRLGQDMEQGVV75HSFFQGRERSRKTIPLVALGLKER



SDS-PAGE Image

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.