

Eukaryotic Plasmin/Antiplasmin Complex (PAP) Instruction Manual

SFPA629Hu61

Homo sapiens (Human)

Source	Eukaryotic expression
Host	293F cell
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted
Predicted Molecular Mass	109.4kDa
Accurate Molecular Mass	120kDa(Analysis of differences refer to the manual)
Residues & Tags	Asn40~Lys491+GGGGS*3+Val98~Arg580 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 5% Trehalose.
Traits	Freeze-dried powder
Purity	> 95%
Isoelectric Point	6.9
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
NQEQVSPLTLLKLGNGQEPGGQTALKSPPGVCSRDPTEQTHRLARAMMAFTADLFSLVAQTSTCPNLILSPLSVALALSHLALGAQNHT
LQRLQQVLHAGSGPCLPHLLSRLCQDLGPGAFRLAARMYLQKGFPIKEDFLEQSEQLFGAKPVSLTGKQEDDLANINQWVKEATEGKIQ
EFLSGLPEDTVLLLLNAIHFQGFWRNKFDPSTLQRDSFHLDEQFTVPVEMMQARTYPLRWFLLEQPEIQVAHFPPKNNMSFVVLVPTHF
EWNVSQVLANLSWDTLHPPLVWERPTKVRLPKLYLKHQMDLVATLSQLGLQELFQAPDLRGISEQSLVSGVQHQSTLELSEVGVEAAA
ATSIAMSRMSLSSFSVNRPFLLFFIFEDTTGLPLFVGSVRNPNPSAPRELKEQQDSPGNKDFLQSLKGFPRGDKLFGPDLKLVPPMEEDY
PQFGSPK
GGGGSGGGSGGGGS
VYLSECKTNGKNYRGTMSKTKNGITCQKWSSTSPHRPRFSPATHPSEGLEENYCRNPDNDPQGPWCYTDDPEKRYDYCDILECEEEECM
HCSGENYDGKISKTMGLECQAWDSQSPHAHGYPISKFPNKNLKNYCRNPDRELRPWCFTTDPNKRWELCDIPRCTTPPPSSGPTYQC
LKGTGENYRGNVAVTVSGHTCQHWSAQTPTHNRTPENFPCKNLDENYCRNPDGKRAPWCHTTNSQVRWEYCKIPSCDSSPVSTEQLAP
TAPPELTPVVQDCYHGDGQSYRGTSSTTTTGKKCQSWSSMTPHRHQKTPENYPNAGLTMNYCRNPDADKGPWCFTTDPNKRWELCDIPRCTTPPPSSGPTYQC
CSGTEASVVAPPVVLLPDVETPSEEDCMFGNGKGYRKRATTVTGTPCQDWAQAEPHRHSIFTPETNPRAGLEKNYCRNPDGDVGGPW
CYTTNPRKLYDYCDVPQCAAPSFDCGKPVQVEPKKCPGR
```

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.