

Eukaryotic Chromogranin A (CHGA) Instruction Manual

SFPB637Hu61

Homo sapiens (Human)

Source	Eukaryotic expression
Host	293F cell
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Secreted, Cytoplasm
Predicted Molecular Mass	16.6kDa
Accurate Molecular Mass	20&27kDa(Analysis of differences refer to the manual)
Residues & Tags	Gln15~Gly148 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Traits	Freeze-dried powder
Purity	> 80%
Isoelectric Point	4.5
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

QVTALPVNSPMNKGDEVMKCI VEVISDTLSKPSMPVVSQECFETLRGDERILSILRHQNL LKELQDLALQGAKERAHQK KHS GFED
ELSEVLENQSSQ AELKEAVEEPSSKDVM EKREDSKEAEKSGEATDG

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

