

# Recombinant Receptor Tyrosine Protein Kinase erbB-2 (ErbB2) Instruction Manual

**SIPB201Hu02**

**Homo sapiens (Human)**

<b>Source</b>	Prokaryotic expression
<b>Host</b>	E.coli
<b>Endotoxin Level</b>	<1.0EU per 1µg (determined by the LAL method)
<b>Subcellular Location</b>	Membrane
<b>Predicted Molecular Mass</b>	42.2kDa
<b>Accurate Molecular Mass</b>	42kDa(Analysis of differences refer to the manual)
<b>Residues &amp; Tags</b>	Thr23~Gly372 with N-terminal His Tag
<b>Buffer Formulation</b>	20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.
<b>Traits</b>	Freeze-dried powder
<b>Purity</b>	> 90%
<b>Isoelectric Point</b>	6.4
<b>Applications</b>	Positive Control; Immunogen; SDS-PAGE; WB.

## SEQUENCE

```
                TQVCTGTD MKLRLPASPE THLDMLRHLY  
QGCQVVQGNL ELTYLPTNAS LSFLQDIQEV QGYVLIAHNQ VRQVPLQRLR  
IVRGTQLFED NYALAVLDNG DPLNNTTPVT GASPGLREL QLRSLTEILK  
GGVLIQRNPQ LCYQDTILWK DIFHKNNQLA LTLIDTNRSR ACHPCSPMCK  
GSRCWGESSE DCQSLTRTVC AGGCARCKGP LPTDCCHEQC AAGCTGPKHS  
DCLAHLHFNH SGICELHCPA LVTYNTDTFE SMPNPEGRYT FGASCVTACP  
YNYLSTDVGS CTLVCP LHNQ EVTAEDGTQR CEKCSKPCAR VCYGLGMEHL  
REVRVTSAN IQEFAGCKKI FG
```

## USAGE

Reconstitute in ddH<sub>2</sub>O to a concentration of 0.1-0.7 mg/mL. Do not vortex.



## STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

## STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

## Image



SDS-PAGE Image

## **[IMPORTANT NOTE]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.