

# Recombinant Interleukin 33 (IL33) Instruction Manual

## SIPB288Mu01

### Mus musculus (Mouse)

<b>Source</b>	Prokaryotic expression
<b>Host</b>	E.coli
<b>Endotoxin Level</b>	<1.0EU per 1µg (determined by the LAL method)
<b>Subcellular Location</b>	Nucleus, Secreted, Cytoplasm, Chromosome
<b>Predicted Molecular Mass</b>	43.1kDa
<b>Accurate Molecular Mass</b>	44kDa(Analysis of differences refer to the manual)
<b>Residues &amp; Tags</b>	Ser109~Ile266 with N-terminal GST Tag
<b>Buffer Formulation</b>	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
<b>Traits</b>	Freeze-dried powder
<b>Purity</b>	> 95%
<b>Isoelectric Point</b>	5.1
<b>Applications</b>	Positive Control; Immunogen; SDS-PAGE; WB.

### SEQUENCE

```
SI QGTSLLTQSP ASLSTYNDQS VSFVLENGCY VINVDDSGKD  
QEQQDQVLLRY YESPCPASQS GDGVDGKKLM VNMSPDKTD IWLHANDKDY  
SVELQRGDVS PPEQAFFVLH KKSSDFVSFE CKNLPGTYIG VKDNQLALVE  
EKDESCNNIM FKLSKI
```

### USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.  
Do not vortex.

### STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

