

Recombinant Anillin (ANLN) Instruction Manual

SIPJ158Hu01

Homo sapiens (Human)

Source	Prokaryotic expression
Host	E.coli
Endotoxin Level	<1.0EU per 1µg (determined by the LAL method)
Subcellular Location	Nucleus, Cytoplasm
Predicted Molecular Mass	50.1kDa
Accurate Molecular Mass	50&35kDa(Analysis of differences refer to the manual)
Residues & Tags	Cys712~Pro1124 with N-terminal His Tag
Buffer Formulation	PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.
Traits	Freeze-dried powder
Purity	> 80%
Isoelectric Point	8.5
Applications	Positive Control; Immunogen; SDS-PAGE; WB.

SEQUENCE

```
CQVNIKQKM QELNNEINMQ QTVIYQASQA LNCCVDEEHG
KGSLEEAEAE RLLLIATGKR TLLIDELNKL KNEGPQRKNK ASPQSEFMPS
KGSVTLSEIR LPLKADFVCS TVQKPDAANY YYLIILKAGA ENMVATPLAS
TSNSLNGDAL TFTTTFTLQD VSNDFEINIE VYSLVQKKDP SGLDKKKKTS
KSKAITPKRL LTSITTKSNI HSSVMASPGG LSAVRTSNFA LVGSYTL SLS
SVGNTKFVLD KVPFLSSLEG HIYLKIKCQV NSSVEERGFL TIFEDVSGFG
AWHRRWCVLS GNCISYWYYP DDEKRKNPIG RINLANCTSR QIEPANREFC
ARRNTFELIT VRPQREDDRE TLVSQCRDTL CVTKNWSAD TKEERDLWMQ
KLNQVLVDIR LWQPDACYKP IGKP
```

USAGE

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

STORAGE

Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

STABILITY

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

Image




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.