

# Recombinant Brother Of CDO Homolog (BOC) Instruction Manual

## SIPL612Mu01

### Mus musculus (Mouse)

<b>Source</b>	Prokaryotic expression
<b>Host</b>	E.coli
<b>Endotoxin Level</b>	<1.0EU per 1µg (determined by the LAL method)
<b>Subcellular Location</b>	Membrane
<b>Predicted Molecular Mass</b>	25.9kDa
<b>Accurate Molecular Mass</b>	23kDa(Analysis of differences refer to the manual)
<b>Residues &amp; Tags</b>	Lys668~Thr866 with N-terminal His Tag
<b>Buffer Formulation</b>	20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.
<b>Traits</b>	Freeze-dried powder
<b>Purity</b>	> 90%
<b>Isoelectric Point</b>	6.3
<b>Applications</b>	Positive Control; Immunogen; SDS-PAGE; WB.

### SEQUENCE

```
          KGI SYKFRVRALN MLGESEPSAP SRPYVMSGYS  
GRVYERPVAG PYITFTDAVN ETTIMLKWMY IPASNNNTPI HGFYIYYRPT  
DSDNDSYKK DMVEGDRYWH SISHLQPETS YDIKMQCFNE GGESEFSNVM  
ICETKARKFS GQPGRPPPLT LAPPQPPPLE TMRPVGTTGA MVARASDLPY  
LIVGVVLGSI VLIIVT
```

### USAGE

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL.  
Do not vortex.

### STORAGE

