SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Recombinant Endoplasmic Reticulum To Nucleus Signalling 1 (ERN1)

Catalog Number: SIPC960Hu01 **Spec: 50ug**

Product Name	Recombinant Endoplasmic Reticulum To Nucleus Signalling 1 (ERN1)
Catalog Number	SIPC960Hu01
Organism Species	Human
Source	Prokaryotic expression
Host	E.coli
Residues	Phe571~Phe832
Tags	N-terminal His Tag
Subcellular Location	Endoplasmic reticulum lumen
Purity	> 90%
Traits	Freeze-dried powder
Buffer formulation	PBS, pH7.4, containing 5% Trehalose.
Original Concentration	400μg/mL
Applications	Positive Control; Immunogen; SDS-PAGE; WB. (May be suitable for use in other assays to be determined by the end user.)
Predicted isoelectric point	6.6
Predicted Molecular Mass	33.3kDa
Accurate Molecular Mass	36kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8ºC for one month.

Aliquot and store at -80ºC for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined

SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

FCPKDVLGHG AEGTIVYRGM FDNRDVAVKR
ILPECFSFAD REVQLLRESD EHPNVIRYFC TEKDRQFQYI AIELCAATLQ
EYVEQKDFAH LGLEPITLLQ QTTSGLAHLH SLNIVHRDLK PHNILISMPN
AHGKIKAMIS DFGLCKKLAV GRHSFSRRSG VPGTEGWIAP EMLSEDCKEN
PTYTVDIFSA GCVFYYVISE GSHPFGKSLQ RQANILLGAC SLDCLHPEKH
EDVIARELIE KMIAMDPQKR PSAKHVLKHP FF

[IDENTIFICATION]




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.