SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Recombinant ATP Binding Cassette Transporter A3 (ABCA3)

Catalog Number: SIPD998Hu01 **Spec: 10ug**

Product Name	Recombinant ATP Binding Cassette Transporter A3 (ABCA3)
Catalog Number	SIPD998Hu01
Organism Species	Human
Source	Prokaryotic expression
Host	E.coli
Residues	Asp1358~Phe1635
Tags	N-terminal His Tag
Subcellular Location	Membrane
Purity	> 95%
Traits	Freeze-dried powder
Buffer formulation	PBS, pH7.4, containing 0.01% Sarcosyl, 5% Trehalose.
Original Concentration	100μg/mL
Applications	Positive Control; Immunogen; SDS-PAGE; WB. (May be suitable for use in other assays to be determined by the end user.)
Predicted isoelectric point	8.5
Predicted Molecular Mass	34.2kDa
Accurate Molecular Mass	35kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH2O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8ºC for one month.

Aliquot and store at -80ºC for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined


SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

DQD VADERTRILA PSPDSLLHTP LIIKELSKVY EQRVPLLAVD
RLSLAVQKGE CFGLLGFNGA GKTTTFKMLT GEESLTSGDA FVGGHRISSD
VGKVRQRIGY CPQFDALLDH MTGREMLVMY ARLRGIPERH IGACVENTLR
GLLLEPHANK LVRTYSGGNK RKLSTGIALI GEPAVIFLDE PSTGMDPVAR
RLLWDTVARA RESGKAIIIT SHSMEECEAL CTRLAIMVQG QFKCLGSPQH
LKSKFGSGYS LRAKVQSEGQ QEALEEFKAF VDLTF

[IDENTIFICATION]



15% SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.