SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Recombinant Eukaryotic Translation Initiation Factor 2C4 (EIF2C4)

Catalog Number: SIPE993Ra01 **Spec: 100ug**

Product Name	Recombinant Eukaryotic Translation Initiation Factor 2C4 (EIF2C4)
Catalog Number	SIPE993Ra01
Organism Species	Rat
Source	Prokaryotic expression
Host	E.coli
Residues	Thr283~Ala597
Tags	N-terminal His Tag
Subcellular Location	Secreted
Purity	> 90%
Traits	Freeze-dried powder
Buffer formulation	pH7.4, containing 0.01% SKL, 5% Trehalose
Original Concentration	200μg/mL
Applications	Positive Control; Immunogen; SDS-PAGE; WB. (May be suitable for use in other assays to be determined by the end user.)
Predicted isoelectric point	8.7
Predicted Molecular Mass	38.9kDa
Accurate Molecular Mass	39kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8ºC for one month.

Aliquot and store at -80ºC for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined


SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

TSPQTLSNLC LKINAKLGGI NNVLVPHQRP SVFQQPVIFL GADVTHPPAG
DGKKPSIAAV VGSMDGHPSR YCATVRVQTS RQEIAQELLF SQEVVQDLTS
MARELLIQFY KSTRFKPTRI IYYRGGVSEG QMKQVAWPEL IAIRKACISL
EEDYRPGITY IVVQKRHHTR LFCADKTERV GKSGNVPAGT TVDSTITHPS
EFDFYLCSHA GIQGTSRPSH YQVLWDDNCF TADELQLLTY QLCHTYVRCT
RSVSIPAPAY YARLVAFRAR YHLVDKDHDS AEGSHVSGQS NGRDPQALAK
AVQIHHDTQH TMYFA

[IDENTIFICATION]



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.