SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Recombinant F-Box And Leucine Rich Repeat Protein 3 (FBXL3)

Catalog Number: SIPJ391Mu01 **Spec: 100ug**

Product Name	Recombinant F-Box And Leucine Rich Repeat Protein 3 (FBXL3)
Catalog Number	SIPJ391Mu01
Organism Species	Mouse
Source	Prokaryotic expression
Host	E.coli
Residues	Met1~Tyr292
Tags	N-terminal His Tag
Subcellular Location	Nucleus, Cytoplasm
Purity	> 97%
Traits	Freeze-dried powder
Buffer formulation	PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.
Original Concentration	200μg/mL
Applications	Positive Control; Immunogen; SDS-PAGE; WB. (May be suitable for use in other assays to be determined by the end user.)
Predicted isoelectric point	7.2
Predicted Molecular Mass	36.7kDa
Accurate Molecular Mass	37kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8ºC for one month.

Aliquot and store at -80ºC for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined


SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

MKRGGRDSDQ DSAEEGTAEK PKRPRTTQER SQPCDWGNLL QDIVLHVFKY
LPLLDRAHAS QVCRNWNQVF HMPDLWRCFE FELNQPATSY LKATHPELIK
QIIKRHSNHL QYVSFKVDSS KESAEAACDI LSQLVNCSLK TLGLISTARP
SFMDLPKSHF ISALTVVFVN SKSLSSLKID DTPVDDPSLK VLVANNSDTL
KLLKMSSCPH VSPAGILCVA DQCHGLRELA LNYHLLSDEL LLALSSEKHV
RLEHLRIDVV SENPGQTHFH TIQKSSWDAF IKHSPKVNLV MY

[IDENTIFICATION]



15% SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.