SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Recombinant C2 Calcium Dependent Domain Containing Protein 3 (C2CD3)

Catalog Number: SIPR996Mu01 **Spec: 100ug**

Product Name	Recombinant C2 Calcium Dependent Domain Containing Protein 3 (C2CD3)
Catalog Number	SIPR996Mu01
Organism Species	Mouse
Source	Prokaryotic expression
Host	E.coli
Residues	Arg2051~Thr2323
Tags	N-terminal His Tag
Subcellular Location	Plasma, Cytoplasm
Purity	> 97%
Traits	Freeze-dried powder
Buffer formulation	PBS, pH7.4, containing 0.01% SKL, 5% Trehalose .
Original Concentration	200μg/mL
Applications	Positive Control; Immunogen; SDS-PAGE; WB. (May be suitable for use in other assays to be determined by the end user.)
Predicted isoelectric point	5.8
Predicted Molecular Mass	33.7kDa
Accurate Molecular Mass	35kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8ºC for one month.

Aliquot and store at -80ºC for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined


SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

RTLNEITTVT DRTSPWSSFM SDMSEVLSPQ PTEVQREGPS CPPEPFPREE
LKVKSSPQKA VSPQPAQGSP SQSGVCEGGA YKIEVEDLAS AKPQPVPSLT FSEAQEGSDS
VGWRASQINQ VRKPMPEMLA ESEAFSSEFS DSSESFETFP LHLPSQSKRE DYKDSPAVRQ
KQVPTGSEVS TRQTLLLPEP VVVPNFFLPP QQLEASLRMI SHSPGLPPAA TTDQDKSEAT
RGALAQRPCR PRPYSIPPNL PEEETRRIAR IFSSQYSKKT EET

[IDENTIFICATION]



15% SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.